

Fig. 2




Fig. 3 Chip holder for 3in x 3in glass chips compatible with standard microscope stages; includes a micro fan for constant "fresh" air, vial holders and electrodes for sample injection






Fig. 4.

## Inlets




T-Inlet classic

Fig. 5(a)



T-Inlet, modified

Fig. 5(b)



T-Inlet, anti-stream

Fig. 5(c)



Inject-Inlet

Fig. 5(d)

Case for "Sheet 6"

## Separation Channel

Fig. 6(a)

Single channel straight




Fig. 6(b)

Single channel meander




Fig. 6(c)

Single channel meander extra long




Fig. 6(d)

Channel bundle parallel, meander




## Evaporators



Single channel


Fig. 7(a)




Funnel-shape

Fig. 7(b)

## Multi Channel Evaporators



Umbel-Shape Fig. 8(a)



Root-Shape Fig. 8(b)




Fig. 8(c)





Fig. 8(d)

1:1 Splitter, rectangular 3-fold      1:1 splitter, rectangular 4-fold




Fig. 9




Chip design with classic T-inlet and medium length meander single channel  
Separator including multi channel root-shape evaporator; all channel  
dimensions a the same (10 $\mu$ m wide and 0.5 $\mu$ m deep)


Fig. 10.



Chip design with inject-inlet including extra long single meander channel for separation ; funnel-evaporator





Chip design including an anti-stream inlet with different angles ( $\alpha$ ,  $\beta$ ) for sample inlet and sample waste, channel dimensions vary between the different regions; bundle of 11 separation channels meandering parallel; evaporator 4-fold 1:1 splitter




Chip design for a three compound synthesis including three umbel-shape  
evaporators and three inject-inlets

fig. 12



Chip design for Immuno-assays including two inject-inlets on two different layers and following "Bessoth-mixer"(Lit); single wide channel evaporator

Fig. 13



design pop02, created 04-04-2000 @ Nils Goedecke

Channel width 110µm after etching, depth 25 µm over the whole structure

Fig. 14(a)




Fig. 14(b)

design popl2, created 14-11-2011 @ Nils Gaedecke




Fig 15(a)

channel width 40 microns for each design

design pop13a by Nils Goedecke 23 June 2000 IC Department of Chemistry

Channel width after etching 60µm; depth 10 µm








Fig. 5(6)



Design Im 01. S I 5, Sep Ch W 10, EVvap Ch W 10 by Nils Goedecke 05.07.2000

Fig. 16 (a)



This layout includes the anti-stream-inlet and a 2.5m separation channel. Theoretically, a channel of this length 10μm wide and 0.1μm deep if running with a  $\eta \sim 40$  has an efficiency of more than 500000 theoretical plates in 10 min run time.

Design 1m 01. S.I. 5, Sep Ch W. 10, EVvap Ch W.10 by Nils Goedecke 05 07 2000

Fig. 16 (b)





Fig. 17(a)

Design Im 02 S I. 5 Sep Ch.W. 10 EV vap Ch.W. 10 by Nils Goedecke 09 11 2000


Goodecke 493400




Design lim 02 S1 5 Sep (h W 10 EVvap (h W 10 by Nils Goedecke 09 11 2000

Fig. 17(b)


a). Design structure  
using CAD package  
and convert to  
machine format




b). Expose  
photoresist  
using DVL II  
system




c). Develop  
photoresist



d). Etch Metal  
Layer



e). Etch Glass



f). Remove Photoresist  
and Metal Layer  
Thermally bond to  
coverplate





Fig. 18



Velocity differences within the channel ( $60 \times 20 \mu\text{m}$ ) for  $10 \mu\text{m}$  latex beads in a pop02 chip driven through evaporation with and without "air condition"; measurement with 50 beads each; The average velocity with the "air condition" switched on is slightly higher than without it – visible in the left shift of the profile.

Fig. 19